

Porous solids for heat storage applications: In-depth material testing by vapor breakthrough measurements

G. Kalies, T. Hähnel, HTW Dresden, Dresden/Germany J. Möllmer, J. Hofmann, Institut für Nichtklassische Chemie e.V., Leipzig/Germany <u>C. Reichenbach</u>, A. Möller, R. Eschrich,

Quantachrome GmbH & Co. KG, Odelzhausen/Germany

# Experimental



- about 80 g of sample where filled in an adsorber and activated at 350 °C for 6 h in 2 L min<sup>-1</sup> nitrogen flow
- **1. Run: Breakthrough of H<sub>2</sub>O** in N<sub>2</sub> carrier gas (80 % RH at 25 °C, total gas flow 4 L min<sup>-1</sup>)
- **Regeneration:** temperature ramp 1 K min<sup>-1</sup> to 130 °C, dwell time 210 min, 2 L min<sup>-1</sup> N<sub>2</sub>-flow  $\rightarrow$  2. Run



**2. Run:** conditions identical to 1. Run **Temperature measurements** along the adsorber axis

### Results



Adsorption capacity / mmol g <sup>-1</sup>	Zeolite 13X BF	Silica Gel
1. Run	18.9	25.9
2. Run	15.4	25.1

- **After Activation:** Integration of 1. Run-breakthrough curves yields adsorption capacities in perfect agreement to water isotherm data measured in a volumetric manner.
- **After Regeneration:** Silica Gel reaches similar adsorption capacities whereas zeolite 13XBF loses about 20 %.



Breakthrough of Water 80 % RH on Silica Gel ----Water 1. Run — Water 2. Run



# Conclusions

- Breakthrough measurements with water vapor provide useful information relevant for adsorptive heat storage applications, e.g. temperature profiles within the adsorber, adsorption capacities, kinetics of adsorption, etc.
- Aging and regeneration studies can be conducted under technically relevant conditions.
- Breakthrough curves can be converted into **>>heat power curves**<</td>characteristic for an open adsorptive heat storage at certain process conditions.
- Energetic considerations (calorimetry) alone miss the impact of dynamical aspects in adsorptive heat storage applications
- Regeneration studies
- Material testing at various conditions
- Direct material comparison at similar conditions
- Temperature profiles within the fixed bed
- Estimation of heat power
- Information on adsorption dynamics



# www.quantachrome.de

Quantachrome GmbH & Co. KG Rudolf-Diesel-Straße 12 D-85235 Odelzhausen

Tel. +49 8134 9324 0 Fax.+49 8134 9324 25 info@quantachrome.de