Characterization of sorption kinetics of CO<sub>2</sub> from N<sub>2</sub>-rich gas mixtures studied by breakthrough experiments on Zeolites

Andreas Möller

Christian Reichenbach, Robert Eschrich



Characterization of

particles · powders · pores

www.3P-instruments.com info@3P-instruments.com



# Application of porous materials as adsorbents

Fine cleaning of gases (e.g. purification of  $H_2$ , natural gas, bio methane...)

Waste air treatment, respiratory protection, solvent recovery, removal of pollutants...)


Gas separation (e.g. air separation...)

Modern and effective materials should have high sorption capacities, high selectivities, and a good kinetic performance.




➡ For such applications, one must consider sorption properties under process-near conditions.










www.3P-instruments.com info@3P-instruments.com




### **Basics – Flow Plot of a Setup for Dynamic Measurements**



- Flow through the regenerated sample with a predefined gas mixture
- Measurement of data at a specified pressure and gas mixture




volume fraction



time


- Determination of technical usable sorption capacity
- Can be used as benchmark for separation performance of adsorbents



time

- Mass transfer coefficient, axial dispersion, shape of isotherm
- Heat effects, heat dissipation
- The time interval of mass transfer zone has to be minimized

### **Saturated Zone**




time

- Determination of saturation capacity
- By assuming of thermodynamic controlled system  $\rightarrow$  Measurement of isotherms possible



Breakthrough curves of 5% CO<sub>2</sub> in N<sub>2</sub> on Zeolites 13X, 5A, 3A (1 bar, 5 l/min (STP), 20°C)



### Qualitative observation of Mass Transfer Zone:

- Zeolite 3A have a spontaneous breakthrough due to too narrow pores (kinetic-steric exclusion)
- Zeolite 5A exhibits a broad mass transfer zone
  - $\rightarrow$  indicates obviously lower kinetic for 5A as 13X
- Both Zeolites, 5A and 13X have quite unsymmetrical breakthrough curves
  - → indicate a **big influence of temperature effects** and shape of isotherms



# Question:


- Is it possible to get reliable kinetic data from such experiments?
- How is the influence of the isotherm shape, temperature effects, flow rates...?

### Answer:

- Yes, but associated with high effort (model of mass- and energy balances is necessary)
  - $\rightarrow$  Simple comparison of slope can be erroneous
- For quantification of temperature effects also a model calculation must be performed!
  - $\rightarrow$  I.e. in some cases heat effects can control nearly the whole curve



## Input Parameters for calculation of breakthrough curves




Kinetics, cycle times, pressure range...

Red:properties of adsorbent/adsorptive systemBlack:properties of adsorber and adsorber wall




# Kinetic considerations - Mass Transfer coefficient $k_{LDF}$



www.3P-instruments.com info@3P-instruments.com





www.3P-instruments.com

info@3P-instruments.com



$$\frac{\partial \overline{q_i}}{\partial t} = k_{LDF} \left( q_i^* - \overline{q_i} \right) \text{ LDF approach } q_i^* = f \left( \text{Isotherms} \right)$$

1)  $k_{LDF,C} = const.$ 

2) 
$$k_{LDF} = \frac{15 \cdot \varepsilon_p \cdot D}{R^2} \frac{c}{q} = k_{LDF,R}^* \frac{c}{q}$$

3)  

$$k_{LDF} = \frac{\rho_p}{A_{SP}} \frac{15 \cdot D}{R^2 \cdot \mu} \frac{1}{1 + \frac{\rho_p}{\varepsilon_p} \frac{\partial q}{\partial c}} = k_{LDF,M}^* \frac{1}{1 + \frac{\rho_p}{\varepsilon_p} \frac{\partial q}{\partial c}}$$
4)  

$$k_{LDF} \sim \frac{D}{R^2} \frac{\partial \ln c}{\partial \ln q} = k_{LDF,D}^* \frac{\partial \ln c}{\partial \ln q}$$

M.G. Plaza et al., *Ind. Eng. Chem. Res*, 55 (2016)
 N.S. Wilkins, A. Rajendran, *Adsorption 25* (2019)
 Mersmann et al., *Chem. Ing. Tech. 55* (1983)
 Darken, *Trans. AIME 175.1* (1948)

#### Note:


$$k_{LDF,C} \neq k_{LDF,M}^* \neq k_{LDF,R}^* \neq k_{LDF,D}^*$$

### **Comparison of LDF values from literature**

 Possible, but used LDF approach must be taken into account!

4... Diffusion in micropores






### Isotherms as input parameter for calculations

- > Knowledge of isotherms necessary due to sharpening and softening of breakthrough curve
- > Favored isotherm (e.g. Type I) shocks adsorption front and softening desorption




# Breakthrough curves for 5% CO<sub>2</sub> on 13X and 5A Zeolite at different flow rates and 20°C



- > Breakthrough curves shifts to shorter breakthrough times and are slightly steeper for higher flow rates
- > Breakthrough curves on Zeolite 13X with smaller Mass Transfer Zone as for Zeolite 5A






- $\blacktriangleright$  Best fits for constant k<sub>LDF</sub> and Darken approach  $\rightarrow$  curves coincides (both are micropore models)
- $\succ$  Best fits for Rajendran and Mersmann  $\rightarrow$  curves coincides (both are macropore models)



### Influence of temperature effects on transport parameter



Adsorption/ Desorption profile on Zeolite 13X

5% CO<sub>2</sub>, 1 bar, 20°C and 0.417 m s<sup>-1</sup>


Modeling with constant  $k_{\text{LDF}}$ 

 $k_{LDF}$ (non-isotherm)= 1.20 min<sup>-1</sup>  $k_{LDF}$ (isotherm) = 0.25 min<sup>-1</sup>

Ratio: 
$$\frac{k_{LDF}(nonisoth)}{k_{LDF}(isoth)} \approx 5$$

- Modeling of both, adsorption and desorption
   Due to influence of isotherm
- Heat effect must be taken into account!








- > More symmetric breakthrough curves for higher temperature (smaller heat effects due to lower loading)
- Breakthrough curves steeper for higher temperatures (faster mass transfer)



## Comparison of breakthrough curves on 5A with different concentrations







1 bar, 20°C and 1.25 m s<sup>-1</sup> 5% CO<sub>2</sub> on 5A k<sub>LDF,C</sub> = 0.4 min<sup>-1</sup>

15% CO<sub>2</sub> on 5A  $k_{LDF,C}$  = 0.8 min<sup>-1</sup>

5% CO<sub>2</sub> on 5A  $k_{LDF,C} = 0.4 \text{ min}^{-1}$ 15% CO<sub>2</sub> on 5A  $k_{LDF,C} = 0.8 \text{ min}^{-1}$ 


1 bar, 20°C and 0.83 m s<sup>-1</sup>

1 bar, 20°C and 0.42 m s<sup>-1</sup> 5% CO<sub>2</sub> on 5A  $k_{LDF,C}$  = 0.4 min<sup>-1</sup> 15% CO<sub>2</sub> on 5A  $k_{LDF,C}$  = 0.5 min<sup>-1</sup>

Strong influence of concentration on breakthrough curves due to non-linearity of isotherm
 Higher LDF-values for higher concentrations



## Comparison of breakthrough curves on 13X and 5A at different flow rates



15% CO2, 1 bar, 20°C and 1.25 m s<sup>-1</sup>15%Zeolite 5A  $k_{LDF,C} = 0.8 \text{ min}^{-1}$ Zeolite

Zeolite 13X  $k_{LDF,C} = 1.2 \text{ min}^{-1}$ 

15%  $CO_2$ , 1 bar, 20°C and 0.83 m s<sup>-1</sup>15%  $CO_2$ Zeolite 5A  $k_{LDF,C} = 0.8 \text{ min}^{-1}$ Zeolite 5/

Zeolite 13X  $k_{LDF,C} = 1.2 \text{ min}^{-1}$ 

15% CO<sub>2</sub>, 1 bar, 20°C and 0.42 m s<sup>-1</sup> Zeolite 5A k<sub>LDF,C</sub> = 0.42 min<sup>-1</sup>


Zeolite 13X  $k_{LDF,C}$  = 1.2 min<sup>-1</sup>

Values same range, and rather independent from flow rate under this conditions
 External film fluid resistance can be neglected



- Model is necessary to get reliable transport parameter → k<sub>LDF</sub> value, with used LDF-Approach
- Strong non-isothermal effects must be considered for evaluation of kinetics → influence of dissipation of heat
- Measurement of desorption part is helpful  $\rightarrow$  influence of isotherm, heat dissipation etc.
- k<sub>LDF</sub> values depend on concentration, total pressure etc. → measurements under same conditions like technical process for up-scaling
- By assuming macropore diffusion resistance of binder → further investigations with different particle sizes from same supplier necessary





Thank you for your attention!

Please visit our website for further information

www.dynamicsorption.com

or contact our partner in US

www.altamirainstruments.com

