Breakthrough Curves of Propane and CO₂ at Different Relative Humidities

Andreas Möller

R. Eschrich, D. Forberg, K. Eisinger, J. Möllmer



Characterization of

particles · powders · pores







Application of porous materials as adsorbents

- Adsorbents for waste air and indoor air cleaning (VOCs like propane)
- Adsorbents for CO₂ separation from waste gas streams
- \blacktriangleright Adsorbents for direct air capture of CO₂
 - For such applications: Humidity plays an important role

Competitve Sorption of H₂O leads to changes in:

- Sorption equilibrium
- Sorption capacity
- Sorption kinetics



www.bluecher.com



Gutknecht et al., Energy Procedia 146 (2018) 129-134; Climeworks AG



Characterization of

Methods



Experimental Setup mix Sorb 0 Q Open Window Options Offline Login Mass Spectrometer Ascii Recalibration 3P 130086 Sample Info Integration Status Table Chart Instrument 13870.7 s (2017-01-19 12:21:36) Waiting for Stop Criteria... (max. 480 min; min. 420 min; TCD threshold [inactive]; TCD stat. [inactive]; Temp. stat. [10 min, <0.1 K]; Pr Vacuum - Exhaust N2 05 bar tum on 6 mL/min 墨 48.7 °C MixSorb L tum on 6 3.042 vol.% mL/min ser-defined (%RH@250 墨 3896.0 tum of 35.5 °C 1904 mL/min 90 14 55.2 °C 90 °C tum off 25.6 °C 4,992 g/h 密 25.0 °C Recirculator tum on 25.00 N2 G ml Unie 300 00 tum off

Flow through the regenerated sample with a predefined gas mixture

Signal detecting (Temp., WLD, etc.) at the outlet for a specified pressure and gas mixture





Microporous Carbon

Microporous Zeolite 4A

Microporous Zeolite 5A

Amin-functionalized microporous Polymer (Lewatit)







Sequence of three Breakthrough curve measurements of propane and H₂O on AC

Segment 1

Segment 2

Segment 3

volume fraction $y(H_2O) = 0.95 \%$, relative humidity approx. 30 % @ 25 °C volume fraction $y(C_3H_8) = 5.00$ %, volume fraction $y(H_2O) = 0.95$ %, relative humidity approx. 30 % @ 25 °C volume fraction $y(C_3H_8) = 5.00 \%$, volume fraction $y(H_2O) = 2.70 \%$, relative humidity approx. 85 % @ 25 °C







Segment 1 (volume fraction $y(H_2O) = 0.95$ %, relative humidity approx. 30 % @ 25 °C



- Breakthrough curve and temperature curves show regular shape
- No distinct condensation observable
- Loading: *q*(H₂O) = 0.64 mmol g⁻¹
- Very good agreement with water isotherm measured on typical volumetric setup





Segment 2 (volume fraction $y(C_3H_8) = 5.00\%$, volume fraction $y(H_2O) = 0.95\%$, relative humidity approx. 30 % @ 25 °C



- Loading: $q(C_3H_8) = 2.84 \text{ mmol g}^{-1}$, good agreement with loading at dry conditions: 2.88 mmol g $^{-1}$
- H_2O isotherm Loading: $q(H_2O) = 0.64 \text{ mmol g}^{-1}$
- Very good agreement with water isotherm measured on typical volumetric setup
- Heat of adsorption by propane is used to desorb water
- Released water can now adsorb in the next adsorbent layer, which is also an exothermal process
- Slight slope in temperature, before big peak observed caused by propane adsorption





Segment 3 (volume fraction $y(C_3H_8) = 5.00$ %, volume fraction $y(H_2O) = 2.70$ %, relative humidity approx. 85 % @ 25 °C)



- Slight displacement of propane, not equilibrated after 5 h → very slow process due to condensation of water in the pores → No replacement of propane on the surface
- By inverting the experiment (1. breakthrough curve with RH 85 % and 2. adding of propane)
 → a lower propane loading (0.79 mmol g⁻¹) is observed
- When adsorbent is pre-saturated with high RH \rightarrow pore filling \rightarrow limited access to adsorption sites





Breakthrough curve of H₂O on Zeolite 4A and Lewatit

(313 K, relative humidity 50 % in N_2 at 0.1 MPa, total flow: 5 Nl min⁻¹)



- Zeolite 4A shows high H₂O sorption capacity compared to Lewatit (due to lower surface area)
- Large heat effects are observed for sorption on zeolite 4A
- Due to lower surface area of Lewatit → the adsorbed amount of H₂O and resulting heat effects are much lower compared to zeolite 4A



1.0



Breakthrough curve on zeolite 4A, 5A and Lewatit – dry vs. wet conditions

(313 K, 2000 ppm CO_2 balanced by N_2 at 0.1 MPa, total flow: 5 Nl min⁻¹, relative humidity 50 %)





- Zeolite 4A and 5A show drastically reduced CO₂ capacities in presence of humidity
- Lewatit shows higher CO₂ capacity in presence of water due to supported reaction of CO₂ with amino-functionality on the surface of the polymer





- Breakthrough experiments are suitable to characterize samples under dynamic conditions in the presence of humidity
- **Presence of humidity** can have large effects on sorption equilibrium and sorption kinetics
- At humidities > 40 %, the loading capacity of propane or CO₂ on AC and zeolite will be reduced
- pre-loaded H₂O will drastically reduce the capacity of propane and CO₂ and thus, the breakthrough time in dynamic experiments
- Competitive sorption in terms of condensation of H₂O (hydrophobic AC, mesoporous adsorbents) can lead to a long equilibration process in dynamic experiments
 → as a result: no equilibrium can be observed from such measurement
- Presence of H₂O is essential in sorption process during amine-functionalized sorbents
 www.3P-instruments.com info@3P-instruments.com www.research.uni-leipzig.de/inc info@inc.uni-leipzig.de







Characterization of particles · powders · pores

> Please visit our website for further information www.dynamicsorption.com

or contact our partner in Brazil

www.acilweber.com.br



Please visit our research and cooperation partner in Germany www.research.uni-leipzig.de/inc

Coming soon as a Digital Event in 2021, April $12^{th} - 15^{th}$:

Adsorption week including 3P Instruments Workshop and INC Symposium