
DYNAMIC WATER SORPTION - EXPERIMENTS AND LESSONS LEARNED



Gunther M. Munz, Sebastian J. Ernst, Andreas Velte, Eric Laurenz, Gerrit Földner

Fraunhofer Institute for Solar Energy Systems ISE

Leipziger Symposium on dynamic sorption

Leipzig, 17th April 2018


www.ise.fraunhofer.de

AGENDA

- Motivation
- Closed applications / pure gas
 - Large Pressure Jumps (LPJ) MeOH/AC Rubotherm TG
 - Large Temperature Jump (LTJ) in kinetic setup
 - Large Pressure Jump (LPJ) in kinetic setup
 - Frequency Response
- Open applications / gas mixtures
 - Dehumidification of (process) air
- Conclusion

Motivation


- Dynamics depend on different processes
 - Mass transfer
 - Reaction speed
 - Heat transfer



Motivation

- Application in focus:



Adsorption heat pumps and chillers → Dynamics = power (density)



Closed / pure gas

LPJ MeOH/AC Rubotherm TG

- Characterization of Activated Carbon extrudate
- Pretreatment / regeneration of the sample
- Pressure Jump after opening a valve
- Detection of uptake (balance)



Closed / pure gas

LPJ MeOH/AC Rubotherm TG


- Characterization of Activated Carbon extrudates

- Results (fitting to model):

- $D_{\text{eff}} = 2-8 * 10^{-5} \text{ m}^2/\text{s}$



- Challenge:

- Keep the adsorbent in isothermal condition
 - No detection of sample temperature → missing information!



Closed / pure gas

Kinetic setup at Fraunhofer ISE



- Large Pressure Jump (LPJ), Large Temperature Jump (LTJ), inert-LTJ, and Small Pressure Jumps (SPJ) are „Standard“ methods



Personnel Information
Name: [Redacted]
Title: [Redacted]
Department: [Redacted]
Phone: [Redacted]
Email: [Redacted]


Safety Information
Warning: [Redacted]
Precautions: [Redacted]
Emergency: [Redacted]



Closed / pure gas

LTJ - Driving temperature differences


- Here: LTJ data for adsorption and desorption of Silicagel (Ø0.9 mm)
- Simple numerical model: Equilibrium data to calculate T_{eqi}
- Temperature difference for heat transfer ΔT_{htTrn} and overall driving temperature difference $\Delta T_{overall}$
- If $\Delta T_{htTrn} \approx \Delta T_{overall}$: Strong heat transfer limitation
- If $\Delta T_{htTrn} \ll \Delta T_{overall}$: Strong mass transfer limitation



Closed / pure gas

LTJ, evaluation of Heat transfer in inert gas

- Here: LTJ data under N_2 atmosphere
- No sorption during temperature jump
- Identification of heat transfer parameters without influence of sorption effects




Closed / pure gas

LPJ - Driving temperature differences





- Here: LPJ data for two samples with a binder based NaY coating (Spray coating)
- Both samples show a strong heat transfer limitation
- But: Adsorption kinetics are very fast
- For the „thick“ sample ($d_{ctg} \approx 450 \mu m$) mass transfer begins to play a role



Closed / pure gas

Identifying heat and mass transfer resistances

- Is it always the heat transfer?
- For different directly crystallized samples (crystallite layer thickness $d_{cryst} = 14 \dots 30 \mu m$) both heat and mass transfer play a role
- Results can help to build optimized heat exchangers



Closed / pure gas

LPJ / LTJ – Modelling: Non-isothermal ads. kinetics (PDE)

Mass balance


$$\frac{\partial c_g}{\partial t} = \frac{\partial}{\partial x} \left(D_{ma} \frac{\partial c_g}{\partial x} \right) + \frac{\dot{n}}{A_{mi}} \frac{4\psi_{mi}}{d_{ma}}$$

$$\frac{\partial c_g}{\partial t} = \frac{1}{r} \frac{\partial}{\partial r} \left(r D_{mi} \frac{c_{ad}}{c_g} \frac{\partial c_g}{\partial r} \right) - \frac{\rho_{crys}^{dry}}{M \psi_{mi}} \frac{\partial X}{\partial t}$$

Heat balance

$$\rho_{comp}^{dry} c_p^{comp} \frac{\partial T_{comp}}{\partial t} = \frac{\partial}{\partial x} \left(\lambda_{comp} \frac{\partial T_{comp}}{\partial x} \right) + h_{ad} \rho_{comp}^{dry} \zeta_{ad} \frac{\partial X}{\partial t}$$

$$h_{ad} = h_{ev} + A$$





Closed / pure gas

LPJ / LTJ – Modelling: Non-isothermal ads. Kinetics (PDE)



Determination of adsorption kinetics parameters

- Experimental values (pressure $p(t)$, temperature $T(t)$, heat flux $Q(t)$) are used to parametrize:
 - Diffusion coefficient
 - thermal conductivity
 - heat transfer coefficients
- Mathematical model of the adsorption process in COMSOL Multiphysics (FEM)
- Diffusion can be described by monodisperse or bidisperse models




Closed / pure gas


Frequency Response Method



- Volume Swing Frequency Response (FR) measurements for in depth investigations of heat and mass as a function of pressure, temperature or loading by model fit [1] in frequency domain.





Closed / pure gas Frequency Response Method



Tab.: Identified coefficients at 11 mbar


$$k_{LDF} = 15 \frac{D}{r^2}$$

Temperature	a (W/m ² K)	D (10 ⁻⁹ m ² /s)
31 °C	230	6
43 °C	280	2
54 °C	220	3
31-54 °C [1]	270	5



Closed / pure gas

Dynamic measurements on different scales





Open applications / gas mixtures

Case Studies

- Dehumidification of (process) air

Open applications / gas mixtures


Quantachrome DynaSorb BT



Open applications / gas mixtures

Case Studies



- Dehumidification of (process) air



Open applications / gas mixtures

Case Studies



■ Dehumidification of (process) air



Open applications / gas mixtures




Case Studies

■ Dehumidification of (process) air






Open applications / gas mixtures

H₂O @ TiAPSO





Open applications / gas mixtures

H₂O @ Al-Fumarate



Dynamic Sorption Measurements



Adsorption – Desorption Cycle



- \dot{V}_{ges} 4500 mL min⁻¹
- 50 % rel. humidity
- $T_{ads} = 32^{\circ}\text{C}$
- $T_{des} = 80^{\circ}\text{C}$

Dynamic Sorption Measurements



Adsorption – Desorption Cycle



- \dot{V}_{ges} 4500 mL min⁻¹
- 50 % rel. humidity
- $T_{ads} = 32^\circ\text{C}$
- $T_{des} = 80^\circ\text{C}$

Open applications / gas mixtures

Adsorption – Desorption Cycle



Conclusion

Closed application:

- Non-isothermal effects are severely influencing derived diffusion parameters
- LPJ and LTJ (+ Inert) are serving well for detecting limitations
- Modelling allows extrapolation to different composite/material and application conditions

Open application:

- Breakthrough experiments deliver the right choice of materials for single and multicomponent adsorption

Thank you for your attention!



Fraunhofer Institute for Solar Energy Systems ISE

Gunther M. Munz, Sebastian J. Ernst

www.ise.fraunhofer.de

gunther.munz@ise.fraunhofer.de; sebastian-johannes.ernst@ise.fraunhofer.de