

# Investigation of Mixed Gas Sorption in Lab-Scale – Experiment and Evaluation

Dr. Andreas Möller

Christian Reichenbach, Robert Eschrich

3P INSTRUMENTS GmbH & Co. KG

2018-06-14

Adsorption Summer School

[www.dynamicsorption.com](http://www.dynamicsorption.com)

## Application of Porous Materials as Adsorbents

Fine cleaning of Gases (i.e. purification of  $H_2$ , natural gas, bio methane...)

Waste air treatment, respiratory protection, solvent recovery, removal of pollutants...)

Gas separation (i.e. Air separation...)

Modern and effective materials should have **high sorption capacities, high selectivities, and a good kinetic performance.**

**For such applications, one must consider gas mixtures and their sorption properties in any case.**



Number of Samples

Application Progress

Synthesis and First Characterization

Determination of Thermodynamic Data

Basic Process Design, Granulation of Adsorbents

Detailed Process Design, Application



Chemists

- BET
- Pore Volume
- Pore Size Distribution



Chemists, Physicists

- Isotherms
- Heat of Adsorption



Chemical Engineers

- Techn. Useable Sorption Capacity
- Gas Mixtures
- Selectivities
- Kinetics for Process Conditions
- Cycle Stability



Engineers

- Process Optimization
- Production

Bench scale, Pilot plants, Industrial Plant

Number of Samples

Application Progress

Synthesis and First Characterization

Determination of Thermodynamic Data

Basic Process Design, Granulation of Adsorbents

Detailed Process Design, Application



Chemists

- BET
- Pore Volume
- Pore Size Distribution



Chemists, Physicists

- Isotherms
- Heat of Adsorption



Chemical Engineers

- Techn. Useable Sorption Capacity
- Gas Mixtures
- Selectivities
- Kinetics for Process Conditions
- Cycle Stability





Engineers

- Process Optimization
- Production

## Textural Properties of Adsorbents:

- BET-Surface
- Pore Size Distribution
- Micropore Volume



**Textural properties allow only limited qualitative statements regarding:**



- **expected sorption capacity (Micropore Volume)**
- **rough assessment of general sorption properties from Pore Size Distribution**

**Textural properties do not allow quantitative statements regarding:**



- **sorption affinity**
- **selectivity**
- **no sufficient information of surface chemistry**

Often different definition of selectivity (separation factor)!

Thermodynamic selectivity and limit of selectivity ( $p \rightarrow 0$ )  $\alpha$

Difference of shape of isotherms and of loadings important [1]

$$\alpha_{CO_2,CH_4} = \frac{Y_{CH_4}}{Y_{CO_2}} \frac{X_{CO_2}}{X_{CH_4}}$$

$$\alpha_{CO_2,CH_4}(p \rightarrow 0) = \frac{H_{CO_2}}{H_{CH_4}}$$

Y... mole ratio in gas phase

X... mole ratio in adsorbed phase

H... Henry constants

Sorbent selection parameter of Rege and Yang [2]

$$S = \frac{\Delta q_{CO_2}}{\Delta q_{CH_4}} \alpha_{CO_2,CH_4}$$

$\Delta q$ ... difference in loading between adsorption and desorption

Kinetic separation faktor  $\beta$

Difference of sorption rates important [3]

$$\beta_{CO_2,CH_4} = \frac{H_{CO_2}}{H_{CH_4}} \cdot \sqrt{\frac{D_{CO_2}}{D_{CH_4}}}$$

H... Henry constants


D... Diffusion coefficients

[1] R.T. Yang, *Gas separation by adsorption processes*, Imperial College Press, London, **1987**

[2] S.U. Rege, R.T. Yang, *Sep. Sci. & Technol.*, **2011**, 36, 3355-3365

[3] D.M. Ruthven, S. Farooq, K.S. Knaebel, *Pressure Swing Adsorption*, Wiley-Verlag, New York **1994**.

## Dependence of partial and total adsorption amounts



General:

$$n_{CO_2,CH_4,total} = FKT(Y_{CO_2}, Y_{CH_4}, p)$$

Investigation along:

THE RED LINE – Case A


$$n_{CO_2,CH_4,total}(p = const.) = FKT(Y_{CO_2}, Y_{CH_4})$$

THE BLUE LINE – Case B

$$n_{CO_2,CH_4,total}(Y_{CO_2}, Y_{CH_4} = const.) = FKT(p)$$


Typical presentation of sorption capacities for binary mixtures

Case A – variable gas composition



$p = \text{const. (1 bar)}$ ,  $\text{CO}_2, \text{CH}_4$  on D55-1.5

Case B – variable pressure



$Y = \text{const. (50:50)}$ ,  $\text{CO}_2, \text{CH}_4$  on D55-1.5




**X-Y-Plot with statement to the composition of adsorbed phase at constant pressure**



**N-p-Plot with statement to the composition of amount of adsorbed at constant composition of the gas phase**


## Typical presentation of sorption capacities for binary mixtures

Case A – variable gas composition



p=const. (1 bar), CO<sub>2</sub>,CH<sub>4</sub> on D55-1.5

Case B – variable pressure



Y=const. (50:50), CO<sub>2</sub>,CH<sub>4</sub> on D55-1.5

➔ **N-Y-Plot with statement to the partial loading at constant pressure**


➔ **N-p-Plot with statement to the composition of amount of adsorbed at constant composition of the gas phase**

## Relationship between loading – mole fraction – selectivity

1. Calculation of **mole fractions of adsorbed phase** from partial loading
2. Calculation of **mole fraction of gas phase** from partial pressures
3. Calculation of **selectivity**
4. Check **plausibility with help of limit for selectivity** (for IAST-Calculations)


$$X_{CO_2} = \frac{n_{CO_2}}{n_{CO_2} + n_{CH_4}} \quad Y_{CO_2} = \frac{p_{CO_2}}{p_{CO_2} + p_{CH_4}}$$

$$\alpha_{CO_2,CH_4} = \frac{Y_{CH_4}}{Y_{CO_2}} \frac{X_{CO_2}}{X_{CH_4}} \quad \alpha_{CO_2,CH_4}(p \rightarrow 0) = \frac{H_{CO_2}}{H_{CH_4}}$$



50% CO<sub>2</sub>, 50% CH<sub>4</sub> on NaMSX, IAST with Toth

- ➔ **Limit of selectivity** can be used to check the results of IAST-Calculations (only for models with Henry range)
- ➔ Often **Limit of selectivity** do not reflect the **selectivity** for the real separation process, therefore a **single consideration is not enough**



## Requirements:

- Knowledge of pure component isotherms

- Multi Component-Langmuir (MCLAI)
- Multi Component-Sips (MCSIPS)
- Multi Component-DSLAI (MCDSLAI)

$$q_i = q_{m,i} \frac{(b_i p_i)^{x_i}}{1 + \sum (b_j p_j)^{x_j}} \quad \text{MCSIPS}$$

- IAST with Langmuir
- IAST with Toth
- IAST with DSLAI

$$\frac{A \cdot \pi}{RT} = \int_0^{p_{i0}} \frac{n}{p} dp = \text{const.}$$


$$n = \text{Isotherm}(p_{i0})$$

- VS-Modell with Wilson

[1] R.T. Yang, *Gas Separation by Adsorption Processes*, Imperial College Press, **1987**

[2] A.L. Myers, J.M. Prausnitz, *AIChE Journal*, **1965**, *11*, 129

[3] T.W. Cochran, R.L. Kabel, R.P. Danner, *AIChE Journal*, *31*, **1985**, *12*, 2075



## Nonlinear Equation System can be solved by Newton Algorithm

Formulation of equation system

$$\begin{aligned} \sum X_i &= 1 \\ \sum Y_i &= 1 \\ Y_i \cdot p &= X_i \cdot p_i^0 \\ \int_0^{p_i^0} \frac{n_i(p)}{p} dp - \frac{\pi A}{RT} &= 0 \end{aligned}$$

Build up of Jacobian matrix  $J$

$$\begin{pmatrix} \frac{\partial f_1}{\partial p_1^0} & \frac{\partial f_1}{\partial p_2^0} & \cdot & \frac{\partial f_1}{\partial \pi} \\ \cdot & \cdot & \cdot & \cdot \\ \cdot & \cdot & \cdot & \cdot \\ \frac{\partial f_N}{\partial p_1^0} & \cdot & \cdot & \frac{\partial f_N}{\partial \pi} \end{pmatrix}$$

Define starting Values (as vector)

$$\begin{pmatrix} p_1^0 \\ \cdot \\ p_i^0 \\ \pi \end{pmatrix}$$

Calculate of negative Functions (as vector)

$$\begin{pmatrix} -FKT_1(p_i^0, \pi) \\ \cdot \\ \cdot \\ -FKT_i(p_i^0, \pi) \end{pmatrix}$$

•

=

Solving of lin. Equation system by **Gaussian elimination**

$$J(p_i^0, \pi) \cdot \Delta(p_i^0, \pi) = -FKTN(p_i^0, \pi)$$

Calculate a new vector

$$p_{i,N+1}^0, \pi_{N+1} = p_{i,N}^0, \pi_N + \Delta(p_i^0, \pi)$$

The whole procedure will be repeated till  $|\Delta(p_{i,N}^0, \pi_N)| < \varepsilon$

**Most Important Problem:**  
**Finding of suitable starting values for convergence!**

**Hint:**

Testing the solution by inserting into equation system

## Example with Langmuir (3 Components)

Formulation of equation system and elimination of  $X_i$

$$\frac{Y_1 \cdot p}{p_1^0} + \frac{Y_2 \cdot p}{p_2^0} + \frac{Y_3 \cdot p}{p_3^0} - 1 = 0$$

$$\int \frac{n_{m1} \cdot b_1}{1 + b_1 \cdot p_1^0} - \Pi = 0$$

$$\int \frac{n_{m2} \cdot b_2}{1 + b_2 \cdot p_2^0} - \Pi = 0$$

$$\int \frac{n_{m3} \cdot b_3}{1 + b_3 \cdot p_3^0} - \Pi = 0$$

Build up of Jacobian matrix  $J$

$$\begin{pmatrix} \frac{Y_1 \cdot p}{(p_1^0)^2} & \frac{Y_2 \cdot p}{(p_2^0)^2} & \frac{Y_3 \cdot p}{(p_3^0)^2} & 0 \\ \frac{n_{m1} \cdot b_1}{1 + b_1 \cdot p_1^0} & 0 & 0 & -1 \\ 0 & \frac{n_{m2} \cdot b_2}{1 + b_2 \cdot p_2^0} & 0 & -1 \\ 0 & 0 & \frac{n_{m3} \cdot b_3}{1 + b_3 \cdot p_3^0} & -1 \end{pmatrix}$$

Define starting Values (as vector)

$$\begin{pmatrix} p_1^0 \\ p_2^0 \\ p_3^0 \\ \pi \end{pmatrix}$$

Calculate of negative Functions (as vector)

$$= -1 \cdot \begin{pmatrix} \frac{Y_1 \cdot p}{p_1^0} + \frac{Y_2 \cdot p}{p_1^0} + \frac{Y_3 \cdot p}{p_1^0} - 1 \\ n_{m1} \cdot \ln(1 + b_1 \cdot p_1^0) - \Pi \\ n_{m2} \cdot \ln(1 + b_2 \cdot p_2^0) - \Pi \\ n_{m3} \cdot \ln(1 + b_3 \cdot p_3^0) - \Pi \end{pmatrix}$$

Solving of lin. Equation system by **Gaussian elimination**

$$J(p_i^0, \pi) \cdot \Delta(p_i^0, \pi) = -FKTN(p_i^0, \pi)$$

$$p_{i,N+1}^0, \pi_{N+1} = p_{i,N}^0, \pi_N + \Delta(p_i^0, \pi)$$

The whole procedure will be repeated till  $|\Delta(p_{i,N}^0, \pi_N)| < \varepsilon$


**Hint:**

In case of Langmuir isotherm an analytical solution exists:

$$p \cdot \left[ \frac{Y_1}{\exp\left(\frac{\Pi}{n_{m1}}\right) - 1} + \frac{Y_2}{\exp\left(\frac{\Pi}{n_{m2}}\right) - 1} + \frac{Y_3}{\exp\left(\frac{\Pi}{n_{m3}}\right) - 1} \right] - 1 = 0$$

## Calculations of mixture data with 3P-Sim – Recommendations for pure components

1. Fitting of pure component data at **same temperature for all components**
2. All data as table - **pressure / bar** (mbar) and **adsorbed amount / mmol g<sup>-1</sup>**
3. **All components must be fitted with same isotherm model**




Typ I: Langmuir, SIPS, Toth, DSLangmuir, DSLangmuirSIPS

Typ II: (Freundlich)


Typ IV, V: DSLangmuirSIPS, (DSLangmuir), (SIPS)

## Calculations of mixture data with 3P-Sim – fitting of pure components



Example: CO<sub>2</sub>, CH<sub>4</sub> on NaMSX at 40°C

## Calculations of mixture data with 3P-Sim – calculation of mixture data



isotherms

➔ **The selected mixture model must include the pure component isotherm model!**  
i.e.: with SIPS only the Multicomponent-SIPS approach is available

## Static-volumetric Method

- Sorption in closed chamber
- Measurement of pressure drop
- Pure component isotherms, **mixed gas sorption hard to realize**
  - Standard characterization (BET, Pore volume, pure component isotherms)




iSorb HP 1

## Dynamic Method

- Sorption in an open system at constant pressure
- Time resolved measurement of effluent gas composition
- **Mainly for mixed gas sorption**, pure component data only with some approaches possible
- Technical processes in „Labscale“
  - Investigations to technical relevant parameters (**selectivity**, sorption kinetic, regenerability, cycle stability ...)



mixSorb L



R.Bazan et al., *Adsorption*, 2011, 17(2), 371-383

- Dosing and **mixing of a predefined gas mixture** in  $V_1$  by circulation pump
- Switch to **adsorption chamber  $V_2$  and continuous mixing** by circulation pump
- Recording of pressure drop and **analysis of gas phase** after achievement of equilibrium
- Calculation of partial loading based on **mass balance**

➔ **Precalculations necessary for a desired equilibrium point**

## Headspace Method (HS):

- **Special case** of a static-volumetric setup
- **Higher sample throughput** as classical volumetry
- Primarily for vapor mixtures suitable

## Principle:

- Regeneration under  $N_2$  at  $150\text{ }^\circ\text{C}$  (AC)
- Sometimes preloading with component 1 (i.e.  $H_2O$ )
- Loading of HS-Vial with sample and **VOC (with  $\mu\text{l}$ -Syringe)**
- **Equilibration for several days** in tempered bath
- **Analysis of gas phase** via GC and performing of a mass balance

➔ **Precalculations necessary for a desired equilibrium point**

M.J.G.Linders et al., *AIChE Journal*, 2001, 47(8), 1885-1892






Headspace GC at INC e.V.

Vials



## Flow plot of a setup for the dynamic method




- Flow through the regenerated sample with a predefined gas mixture
- Measurement of data at a specified pressure and gas composition → **Advantage over static volumetry, no precalculations necessary**

A.Möller et al., *Adsorption*, 2016, 23(2-3), 197-209

Static-Volumetric	Headspace Method	Dynamic Method
➤ Closed Chamber	➤ Closed Chamber	➤ Open System
➤ no carrier gas necessary	➤ no carrier gas necessary	➤ carrier gas necessary, depending on the routine
➤ Suitable for gas mixtures and vapors	➤ Mainly for vapors	➤ Suitable for gas mixtures and vapors
➤ Experimentally complex	➤ <b>Precalculations necessary</b>	➤ <b>No precalculations</b> necessary
➤ <b>Precalculations necessary</b>	➤ Sequenced experiments hard to perform	➤ Sequenced experiments easy to perform
➤ Sequenced experiments hard to perform		➤ <b>Experiments along the „Lines“ possible</b>

Static-Volumetric	Headspace Method	Dynamic Method
➤ Closed Chamber	➤ Closed Chamber	➤ Open System
➤ no carrier gas necessary	➤ no carrier gas necessary	➤ carrier gas necessary, depending on the routine
➤ Suitable for gas mixtures and vapors	➤ Mainly for vapors	➤ Suitable for gas mixtures and vapors
➤ Experimentally complex	➤ <b>Precalculations necessary</b>	➤ <b>No precalculations necessary</b>
➤ <b>Precalculations necessary</b>	➤ Sequenced experiments hard to perform	➤ Sequenced experiments easy to perform
➤ Sequenced experiments hard to perform		➤ <b>Experiments along the „Lines“ possible</b>


## Balancing of breakthrough experiments



- Widely used approach and simplification for high diluted systems:

$$\dot{V}_{\text{out}}(t) = \dot{V}_{\text{in}}(t) = \text{const.}$$

- ➔ **Small changes in concentration minimize this error**
- ➔ **High diluted systems easier to evaluate**



- **binary mixture:**

- CO<sub>2</sub>/He (non-adsorbable carrier gas)


- → Pure component equilibria

- CO<sub>2</sub>/CH<sub>4</sub> (adsorbable carrier gas)

- → Preloading of sample with pure CH<sub>4</sub>

- → Incomplete Determination of the system  
(evaluation mostly simple)

- → Partial loading for CO<sub>2</sub> (mixture sorption data)



- **ternary mixture:**

- CO<sub>2</sub>/CH<sub>4</sub>/He (non-adsorbable carrier gas)

- Displacement of less adsorbed component

- → Partial Desorption, Role-Up Effects

- → Complete Determination (evaluation complex)


- CO<sub>2</sub>/CH<sub>4</sub>/N<sub>2</sub> (adsorbable carrier gas)

- → Preloading of sample with pure N<sub>2</sub>


- → Partial loading for CO<sub>2</sub> and CH<sub>4</sub>

- → Incomplete ternary mixture data


- Calculation of mixture data desired (Y/N)?
  - Pure component isotherms necessary
- Is a **complete determination of the system** desired?
  - **determination of all partial loadings**, diluting with Helium-carrier gas (Y/N)
- Definition of total flow, concentration, measurement temperature etc.
  - Sample must be under „thermodynamic control“ (always)
- Depending on **concentration range** one should consider:
  - Calibration of suitable analytic technique (always)
  - Check if any approach is valid, i.e. constant gas velocity ...(always)
- Sample preparation and definition of preparation conditions
  - Temperature, carrier gas (always)
- Build up of a measurement routine
  - pressurization, Helium or Adsorptive 1 (Helium for complete determination)
- Evaluation of the experiment




- Calculation of mixture data desired (Y/N)?
  - Pure component isotherms necessary
- Is a **complete determination of the system** desired?
  - **determination of all partial loadings**, diluting with Helium-carrier gas (Y/N)
- Definition of total flow, concentration, measurement temperature etc.
  - **Sample must be under „thermodynamic control“** (always)
- Depending on **concentration range** one should consider:
  - Calibration of suitable analytic technique (always)
  - Check if any approach is valid, i.e. constant gas velocity ...(always)
- Sample preparation and definition of preparation conditions
  - Temperature, carrier gas (always)
- Build up of a measurement routine
  - pressurization, Helium or Adsorptive 1 (Helium for complete determination)
- Evaluation of the experiment



- Calculation of mixture data desired (Y/N)?
  - Pure component isotherms necessary
- Is a **complete determination of the system** desired?
  - **determination of all partial loadings**, diluting with Helium-carrier gas (Y/N)
- Definition of total flow, concentration, measurement temperature etc.
  - **Sample must be under „thermodynamic control“** (always)
- Depending on **concentration range** one should consider:
  - Calibration of suitable analytic technique (always)
  - Check if any approach is valid, i.e. constant gas velocity ...(always)
- Sample preparation and definition of preparation conditions
  - Temperature, carrier gas (always)
- Build up of a measurement routine
  - pressurization, Helium or Adsorptive 1 (Helium for complete determination)
- Evaluation of the experiment




- Calculation of mixture data desired (Y/N)?
  - Pure component isotherms necessary
- Is a **complete determination of the system** desired?
  - **determination of all partial loadings**, diluting with Helium-carrier gas (Y/N)
- Definition of total flow, concentration, measurement temperature etc.
  - **Sample must be under „thermodynamic control“** (always)
- Depending on **concentration range** one should consider:
  - Calibration of suitable analytic technique (always)
  - Check if any approach is valid, i.e. constant gas velocity ...(always)
- Sample preparation and definition of preparation conditions
  - Temperature, carrier gas (always)
- Build up of a measurement routine
  - pressurization, Helium or Adsorptive 1 (**Helium for complete determination**)
- Evaluation of the experiment



## Task: Investigation of a binary system


Activated Carbon, CO<sub>2</sub> (25%), CH<sub>4</sub> (75%), complete determination at 5 bar

1. Weighting the sample and sample preparation at 120°C, He-flow 200 ml min<sup>-1</sup> (STP)
2. Definition of partial pressures: 1,25 bar CO<sub>2</sub>; 3,75 bar CH<sub>4</sub>; 5 bar He; Σ10 bar
3. Gas flows: 1l min<sup>-1</sup> (STP) He, 0,75 l min<sup>-1</sup> (STP) CH<sub>4</sub>, 0,25 l min<sup>-1</sup> (STP) CO<sub>2</sub>
4. Pressurization with Helium up to 10 bar
5. Start of measurement by simultaneous dosing of CO<sub>2</sub> and CH<sub>4</sub> in Helium
6. Recording of effluent gas composition via MS (all components!)
7. After breakthrough, regeneration of sample for determination of activated mass



## Task: Investigation of a binary system

Activated Carbon, CO<sub>2</sub> (25%), CH<sub>4</sub> (75%), complete determination at 5 bar



### Result:



Breakthrough curve with “Role-Up” Effect

- Includes all partial loadings
- Reference on  $p_{\text{CO}_2}=0,25$  bar and  $P_{\text{CH}_4}=3,75$  bar
- Mole fraction:  $y_{\text{CO}_2}=0,25$ ;  $y_{\text{CH}_4}=0,75$
- Helium will not be considered!

### Integration of areas:

- $n(\text{CO}_2)$ ;  $n(\text{CH}_4)$ ;  $n(\text{total})$ ;  $\alpha$

5% CO<sub>2</sub> 15% CH<sub>4</sub> in He at 20°C, 5 bar, 2500 ml min<sup>-1</sup> (STP) on D 55/1.5



Thermodyn. Model  
(Multi component-SIPS)


$$Y_{\text{CO}_2} = 0,25$$

$$Y_{\text{CH}_4} = 0,75$$

$$P_{\text{total}} = 1 \text{ bar}$$

$$\alpha_{\text{Model}} = 3,53$$

A complete determined System!



Dynamic experiment  
(determination of all partial loadings)



$$n_{\text{CO}_2} = 0,86 \text{ mmol g}^{-1}$$

$$n_{\text{CH}_4} = 0,71 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Experiment}} = 3,63$$

→ System is complete determined, good agreement confirmed ideal behavior.

5% CO<sub>2</sub> 95% N<sub>2</sub> at 40°C, 5 bar, 2000 ml min<sup>-1</sup> on D 55/1.5



Thermodyn. Model  
(Multi component-SIPS)


$$Y_{\text{CO}_2} = 0,05$$

$$Y_{\text{N}_2} = 0,95$$

$$P_{\text{total}} = 5 \text{ bar}$$

$$\alpha_{\text{Model}} = 15,7$$

System is not complete determined!



Dynamic Experiment  
(Determination of  
on partial loading)


$$n_{\text{CO}_2} = 0,57 \text{ mmol g}^{-1}$$

$n_{\text{N}_2}$  = not determined

$\alpha_{\text{Experiment}}$  = not determined

→ System is only incomplete determined. Model is confirmed by experiment.

## Sequence of several breakthrough curves on Activated Carbon D 55-1.5



### Conditions:


- 20 °C, 2 L min<sup>-1</sup>
- 10 bar (pressurization with N<sub>2</sub>)
- Concentrations: 5 % CO<sub>2</sub> - 80 % CO<sub>2</sub> in N<sub>2</sub>

### Procedure:

- Start further breakthroughs after equilibrium before
- **Integration** and **summation** results in partial loading data of CO<sub>2</sub>
- Volume ratio and total pressure defines the partial pressure of CO<sub>2</sub>

→ **Mixed isotherm data of CO<sub>2</sub> in N<sub>2</sub>**

## Measured partial loading data for CO<sub>2</sub> on Activated Carbon D 55-1.5 at 10 bar



- Dynamic measured data (red)
- **IAST**-Model (Ideal Adsorbed Solution Theory) based on pure component isotherms (lines)
- Mixture of CO<sub>2</sub> and N<sub>2</sub> shows **ideal** behavior on AC




Determination of partial loading data of CO<sub>2</sub> on AC D 55-1.5 by performing sequentially experiments along constant total pressure.

$$n_{\text{CO}_2, \text{N}_2, \text{total}}(p = \text{const.}) = FKT(Y_{\text{CO}_2}, Y_{\text{N}_2})$$



## Sequence of several breakthrough curves on Activated Carbon D 55-1.5

- Assumption of non-adsorbable Gas, i.e. **Helium** as carrier gas will lead to pure component equilibria.
- Other, **adsorbable gases** results in mixture sorption data.





$\text{CO}_2/\text{He}$

$\text{CO}_2/\text{N}_2$



1% EtOH 2,5% H<sub>2</sub>O (approx. 80% RH) in N<sub>2</sub> at 25°C, 1 bar, 4000 ml/min on NaMSX




Thermodyn. Model  
(IAST-DSLAI)

$$Y_{\text{EtOH}} = 0,322$$

$$Y_{\text{H}_2\text{O}} = 0,678$$

$$P_{\text{total}} = 0.033 \text{ bar}$$

$$\alpha_{\text{Model}} > 20000$$



Dynamic experiment  
(determination of all  
partial loadings)



$$n_{\text{EtOH}} = \text{ca. } 0 \text{ mmol g}^{-1}$$

$$n_{\text{H}_2\text{O}} = 12,9 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Experiment}} = \text{not determined } (>20)$$

➔ **Very high selectivity of model confirmed by experiment**

1% EtOH 2,5% H<sub>2</sub>O (approx. 80% RH) in N<sub>2</sub> at 25°C, 1 bar, 4000 ml/min on D 55/1.5




Thermodyn. Model  
(IAST-DSLAIIPS)

$$Y_{\text{EtOH}} = 0,286$$

$$Y_{\text{H}_2\text{O}} = 0,714$$

$$P_{\text{total}} = 0.035 \text{ bar}$$

$$\alpha_{\text{Model}} = 2,9$$



Dynamic experiment  
(determination of **all partial loadings**)


$$n_{\text{EtOH}} = 3,82 \text{ mmol g}^{-1}$$

$$n_{\text{H}_2\text{O}} = 3,11 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Experiment}} = 3,1$$

➔ **Deviation between model and experiment** ➔ no prediction, experiment necessary!

15% Propane 45% Propylene in He at 25°C, 5 bar, 1000 ml/min on A1, A2, A3



Sample: A1


$$Y_{\text{Propane}} = 0,25$$

$$Y_{\text{Propylene}} = 0,75$$

$$n_{\text{Propane}} = 0,03 \text{ mmol g}^{-1}$$

$$n_{\text{Propylene}} = 0,06 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Propylene}} = 0,67$$



Sample: A2


$$Y_{\text{Propane}} = 0,25$$

$$Y_{\text{Propylene}} = 0,75$$

$$n_{\text{Propane}} = 0,2 \text{ mmol g}^{-1}$$

$$n_{\text{Propylene}} = 1,7 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Propylene}} = 2,83$$



Sample: A3

$$Y_{\text{Propane}} = 0,25$$

$$Y_{\text{Propylene}} = 0,75$$

$$n_{\text{Propane}} = \text{ca. } 0 \text{ mmol g}^{-1}$$

$$n_{\text{Propylene}} = 2,1 \text{ mmol g}^{-1}$$

$$\alpha_{\text{Propylene}} = \text{not determined } (>20)$$

Selectivity

- ➔ Statements on selectivity also possible without modeling from pure component data
- ➔ Determination of sorption capacities and selectivities, recording of kinetic

## ➤ Relevance for materials / science and production

- Assessment of porous materials
- R&D of tailor made materials for certain separation issues
- Evaluation of Influence of some secondary components



Zeolite for gas purification  
Source: [www.cwk-bk.de](http://www.cwk-bk.de)

## ➤ Relevance for engineering / Process design

- Modeling and Calculation of fixed bed adsorbers
- Determination of Kinetics (Linear driving force constants)
- Calculation of breakthrough curves, upscaling by modeling
- Calculation of PSA and TSA processes



Process gas purification plant  
source [www.silica.de](http://www.silica.de)

## ➤ Relevance for materials / science and production

- Assessment of porous materials
- R&D of tailor made materials for certain separation issues
- Evaluation of Influence of some secondary components



Zeolite for gas purification  
Source: [www.cwk-bk.de](http://www.cwk-bk.de)


## ➤ Relevance for engineering / Process design

- Modeling and Calculation of fixed bed adsorbers
- Determination of Kinetics (Linear driving force constants)
- Calculation of breakthrough curves, upscaling by modeling
- Calculation of PSA and TSA processes



Process gas purification plant  
source [www.silica.de](http://www.silica.de)

## Input Parameters





Cycle duration, pressure range...

**Red:** properties of adsorbent/adsorptive system  
**Black:** properties of adsorber and adsorber wall

# Modeling of Breakthrough Curves

## Kinetic considerations - Mass Transfer coefficient $k_{LDF}$




### Simplification



## Mass Balance

$$\underbrace{\frac{\partial C_{g,i}}{\partial t} + \frac{1-\varepsilon}{\varepsilon} \rho_{part} \frac{\partial \bar{q}_i}{\partial t}}_{\text{accumulation}} - \underbrace{Dax \frac{\partial^2 C_{g,i}}{\partial z^2}}_{\text{dispersion}} + \underbrace{u \frac{\partial C_{g,i}}{\partial z} + C_{g,i} \frac{\partial u}{\partial z}}_{\text{convection}} = 0$$

$$\frac{\partial \bar{q}_i}{\partial t} = k_{LDF} (q_i^* - \bar{q}_i) \quad \text{LDF approach} \quad q_i^* = f(\text{Isotherms})$$



## Energy Balances

$$\underbrace{\frac{1-\varepsilon}{\varepsilon} \rho_{part} \sum_{i=1}^n \frac{\Delta H_i}{M_i} \frac{\partial \bar{q}_i}{\partial t}}_{\text{generation}} - \underbrace{\lambda \frac{\partial^2 T_g}{\partial z^2}}_{\text{dispersion}} + \underbrace{\left( \frac{1-\varepsilon}{\varepsilon} \rho_{part} cps + \rho_g cpg \right) \frac{\partial T_g}{\partial t}}_{\text{accumulation}} + \underbrace{u \rho_g cpg \frac{\partial T_g}{\partial z} + T_g \rho_g cpg \frac{\partial u}{\partial z}}_{\text{convection}} + \underbrace{4 \frac{h_w}{d_i} (T_g - T_w)}_{\text{transfer to wall}} = 0$$



$$\underbrace{\alpha_w h_w (T_g - T_w)}_{\text{gas to wall}} - \underbrace{\rho_w cpw \frac{\partial T_w}{\partial t}}_{\text{accumulation}} - \underbrace{\alpha_{wL} U_g (T_w - T_{Env})}_{\text{wall to environment}} = 0$$

## Equation for velocity / overall mass balance (isothermal)


$$\frac{\partial u}{\partial z} + \underbrace{\frac{RT}{p} \rho_{part} \frac{1-\varepsilon}{\varepsilon} \sum_{i=1}^n \frac{1}{M_i} \frac{\partial \bar{q}_i}{\partial t}}_{\text{Change by adsorption}} + \underbrace{\frac{1}{p} \frac{\partial p}{\partial t}}_{\text{Change by compression}} = 0$$

\* A.M. Ribeiro et. al, *Chem. Eng. Sci.* 63 (2008)  
 \* D.M. Ruthven, *Principles of Adsorption* (1984)  
 \* M.S. Shafeeyan et. al, *Chem. Eng. Res. Des.* 92 (2014)

## Replacement effects





Input pure component Isotherms



Fitting  $k_{LDF}$

Fitting Heat Transfer Coeff.



5% CO<sub>2</sub> 15% CH<sub>4</sub> in He at 20°C, 5 bar, 2500 ml/min on D 55/1.5

## Replacement effects – Temperature profiles

Use of integral heat of adsorption with Sips model (from isotherms):

$$q_{s,i} = q_{\max,i} \cdot \frac{(K_i \cdot c_i)^i}{1 + \sum_{j=1}^n (K_j \cdot c_j)^j}$$

$$K_i = K_{i,0} \cdot \exp\left(\frac{Q}{R} \left(\frac{1}{T} - \frac{1}{T_0}\right)\right)$$

$$q_{\max,i} = q_{\max,i,0} \cdot \exp\left(\chi_i \left(1 - \frac{T}{T_0}\right)\right)$$

$$t_i = t_{i,0} + \alpha_i \left(1 - \frac{T_0}{T}\right)$$


Input Q in energy balance

Input  $h_w$ ,  $U_g$  in energy balance

- Q (equivalent to  $-\Delta H_{\Theta=0,5}$ )
- $Q_{CO_2}$ : 15.9 kJ mol<sup>-1</sup>
- $Q_{CH_4}$ : 11.9 kJ mol<sup>-1</sup>

### Fitting Heat Transfer Coeff.


- $h_w$  : ~ 30 W m<sup>-2</sup> K<sup>-1</sup> (Gas, Fixed Bed/Wall)
- $U_g$  : ~ 400 W m<sup>-2</sup> K<sup>-1</sup> (Wall/Liquid)



5% CO<sub>2</sub> 15% CH<sub>4</sub> in He at 20°C, 5 bar, 2500 ml/min on D 55/1.5

- CH<sub>4</sub> induced higher temperature effect
- Model can describe temperature profiles qualitatively
- Underestimation of temperature peaks
- Experiment shows mostly sharper temperature profiles
  - differences due to simplification of no radial gradients
  - radial gradients in experiment expected due to external liquid cooling

## Example: CO<sub>2</sub> Adsorption on D55/1.5 in N<sub>2</sub>-Carrier Gas



### Observations:

- 1) Desorption curve flatter than Adsorption curve
- 2) Desorption time higher than Adsorption time
- 3) Adsorption time 5.78 min ( $c_{out} < 0.2\%$ )


### Questions concerning:

- 1) Kinetic parameter ( $k_{LDF}$ )
- 2) Total pressure during each step
- 3) Adsorption/Desorption times
- 4) Purge flow during Desorption

**Adsorption:** 5% CO<sub>2</sub> in N<sub>2</sub> at 40°C, 5 bar, 2000 ml/min on D 55/1.5

**Desorption:** Purging with 2000 ml/min N<sub>2</sub> at 40°C, 5 bar

## Regeneration / PSA



5% CO<sub>2</sub> 95% N<sub>2</sub> at 40°C, 5 bar, 2000 ml/min on D 55/1.5

### Model after Fitting

- Isotherms (MCSIPS)
- Kinetic parameter ( $k_{LDF}$ )
- Heat Transfer Parameter

→ Model can consider slower Desorption due to curved isotherm

### Parameter from Experiment:

- Adsorption time 5.78 min
- Adsorption pressure 5 bar
- Feed Flow 2000 ml/min
- Purge Flow 2000 ml/min pure N<sub>2</sub>

### General Requirements for PSA:

- Purge Flow 500 ml/min pure N<sub>2</sub>
- Desorption in counter current flow
- Max. CO<sub>2</sub> content in product 1%


### Question concerning:

**1) Desorption pressure?**

## Regeneration / PSA

### Cycle times for modeling:

- Adsorption time 5.78 min
- Desorption time 5.03 min
- Calculating 5 Cycles



### Cycle times for Experiment:

- Adsorption time 5.78 min @ 5 bar
- Blow Down time ~ 0.25 min
- Desorption time 4.75 min
- Pressurization to 4.6 bar with N<sub>2</sub>
- Pressurization from 4.6 bar to 5 bar with Feed
- Measurement of 5 cycles

Calculations with  $p_{DES} = 1$  bar


Predictions by modeling:

**Regeneration conditions not strong enough  
→ CO<sub>2</sub> impurity in effluent flow increases  
from cycle to cycle to ~ 3 %**

## Regeneration / PSA

### Cycle times for modeling:

- Adsorption time 5.78 min
- Desorption time 5.03 min
- Calculating 5 Cycles



### Cycle times for Experiment:

- Adsorption time 5.78 min @ 5 bar
- Blow Down time ~ 0.25 min
- Desorption time 4.75 min
- Pressurization to 4.6 bar with N<sub>2</sub>
- Pressurization from 4.6 bar to 5 bar with Feed
- Measurement of 5 cycles

Calculations with  $p_{DES} = 1$  bar

Predictions by modeling:


**Regeneration conditions not strong enough  
→ CO<sub>2</sub> impurity in effluent flow increases  
from cycle to cycle to ~ 3 %**

Predictions were confirmed by experiment

## Regeneration / VPSA

### Cycle times for modeling:

- Adsorption time 5.78 min
- Desorption time 5.03 min
- Calculating 5 Cycles



### Cycle times for Experiment:

- Adsorption time 5.78 min @ 5 bar
- Blow Down time ~ 0.25 min
- Desorption time 4.75 min
- Pressurization to 4.6 bar with N<sub>2</sub>
- Pressurization from 4.6 bar to 5 bar with Feed
- Measurement of 5 cycles

Calculations with  $p_{DES} = 0.5$  bar

Predictions by modeling:

**Regeneration conditions good enough**

**→ CO<sub>2</sub> impurity in effluent flow increases from cycle to cycle, but still below target (<1%)**


## Regeneration / VPSA

### Cycle times for modeling:

- Adsorption time 5.78 min
- Desorption time 5.03 min
- Calculating 5 Cycles

### Cycle times for Experiment:

- Adsorption time 5.78 min @ 5 bar
- Blow Down time ~ 0.25 min
- Desorption time 4.75 min
- Pressurization to 4.6 bar with N<sub>2</sub>
- Pressurization from 4.6 bar to 5 bar with Feed
- Measurement of 5 cycles



Calculations with  $p_{DES} = 0.5$  bar

Predictions by modeling:

**Regeneration conditions good enough**

**→ CO<sub>2</sub> impurity in effluent flow increases from cycle to cycle, but still below target (<1%)**

Predictions were confirmed by experiment


## Regeneration / VPSA

### Cycle times for modeling:

- Adsorption time 5.78 min
- Desorption time 5.03 min
- Calculating 5 Cycles

### Cycle times for Experiment:

- Adsorption time 5.78 min @ 5 bar
- Blow Down time ~ 0.25 min
- Desorption time 4.75 min
- Pressurization to 4.6 bar with N<sub>2</sub>
- Pressurization from 4.6 bar to 5 bar with Feed
- Measurement of 5 cycles



### But: modeling differs from experiment!



- Cycle Steps in modeling strong simplified
- Variations experiment from model mainly in desorption part

- Modeling can help to reduce experimental effort
- final evaluation only by experiment!

## $C_3H_8$ removal from $CH_4$ (partial pressure range 0.01 bar and 0.50 bar)

Selection of activated carbon with different BET-Surfaces, but from same raw material

- AC 1 BET  $\sim 1800 \text{ m}^2/\text{g}$
- AC 2 BET  $\sim 1300 \text{ m}^2/\text{g}$



Often AC with higher BET will be selected by user which is not always the best decision!


- according to isotherms AC 2 is better for low  $C_3H_8$  concentrations
- for high  $C_3H_8$  concentrations AC 1 is better

## $C_3H_8$ removal from $CH_4$ (partial pressure range 0.01 bar and 0.50 bar)


Selection of activated carbon with different BET-Surfaces, but from same raw material

- AC 1 BET  $\sim 1800 \text{ m}^2/\text{g}$
- AC 2 BET  $\sim 1300 \text{ m}^2/\text{g}$

breakthrough of 1%  $C_3H_8$  in  $CH_4$  at  $40^\circ\text{C}$ , 1 bar



breakthrough of 50%  $C_3H_8$  in  $CH_4$  at  $40^\circ\text{C}$ , 1 bar






**Breakthrough experiments and simulations very sensitive for low concentrations!**

→ Observations made from isotherms were confirmed by dynamic experiments and calculations

## Calculation of Constant Pattern Profiles

For favored isotherms (Type I-Isotherms) a Constant Pattern Behavior can occur

- Shape of breakthrough will not change for longer elongation times or adsorber heights, respectively
- Based on compensation of flattening and rising effects
- Height for Constant Pattern =  $f(\text{Shape of Isotherm, Dispersion, Kinetics})$




5% CO<sub>2</sub> in He at 40°C, 5 bar, 1000 ml/min on D 55/1.5

Experiment carried out at 20 cm, simulations were performed for different heights

- Slopes at  $C/C_0=0.5$  were used to evaluate steepness of breakthrough curves
- Constant Pattern Behavior can be expected above 20 cm bed height

- Equilibria of mixtures and selectivities are **important** for applications
  - Knowledge of such data are highly valuable and of great interest
  - however much **more difficult to measure**
    - There is a lack of data in literature
- Mixture equilibria can not be calculated by textural properties but **partial predictive from pure component isotherms** based on classical thermodynamic Models.
- **Predictions have some limits! → Number of components, vapors, non-ideal systems..., experiments are mandatory**
- Mixture Models necessary for calculations of dynamic processes





Thank you for your attention!

Please visit our website for further information

[www.dynamicsorption.com](http://www.dynamicsorption.com)

via E-Mail on

[info@3p-instruments.com](mailto:info@3p-instruments.com)



You are invited!



LEIPZIGER SYMPOSIUM  
on dynamic sorption

2019

»Advanced Sorbent Materials on the Way to Application«

May 14th 2019 in Leipzig (Germany)  
[www.dynamicsorption.com/symposium](http://www.dynamicsorption.com/symposium)

*Dynamic sorption*